

DATASHEET

IPA Spray

FEATURES AND BENEFITS

IPA 30% DI water Filtered to 0.2 micron, White 1 Litre bottles with 28ml trigger 70% spray, double bagged 6 bottles per box, can be gamma irradiated with not less than .25kGy, specific gravity 0.872 @20degC 2 years shelf life from date of manufacture

- 70% IPA 30% D.I Water 0.2mu filtered
- Ideal for the disinfection of product contact surfaces
- Broad Spectrum Activity
- Fast contact time
- Double bagged

APPLICATIONS

Cleanrooms, Hospitals and Laboratories.

FLUID			
IPA	70%		
MATERIAL	DI Water 30%		
	Filtered to 0.2 Micron		
PACKAGING	White 1 litre bottles with 28ml Trigger Spray		
DOUBLE BAGGED	6 bottles per box		
IRRADIATION	Can be Gamma Irradiated with not less than 25kGy		
SPECIFIC GRAVITY	0.871 @ 20degC		
SHELF LIFE	2 years from date of manufacture		

PRODUCT CODE	SIZE		
600-0801	Trigger Spray	Non-Sterile	1 x 1 Litre
600-0802G	Trigger Spray	Sterile	1 x 1 Litre
600-0803	RTU solution	Non-Sterile	1 x 2.5 Litre



AL JAZEIRA ST., AL RIGGA RD. **DUBAI, PO BOX 171315** UNITED ARAB EMIRATES PHONE: +971 4 2959818 FAX: +971 4 2501223

E-MAIL: SALES@WIDACO.COM















